



CLEANING AND DISINFECTION IN HEALTHCARE ENVIRONMENTS

2 in 1 Detergent and Disinfectant

For fixed surfaces
and non-critical items

5th Generation
Polymeric Biguanide and
Quaternary Ammonium



5L Drum Model

Sprayer Model 750mL



READY TO USE



FAST ACTION



BIODEGRADABLE



ECONOMIC

FAST ACTION IN

2 min

- **Superfungus:**
 - Candida auris
- **Gram-negative bacteria:**
 - Escherichia coli
 - Salmonella choleraesuis
 - Pseudomonas aeruginosa
 - Klebsiella pneumoniae carbapenemase (KPC)
 - Acinetobacter baumannii
- **Gram-positive bacteria:**
 - Enterococcus hirae
 - Staphylococcus aureus
 - Staphylococcus aureus MRSA
 - Enterococcus faecalis VRE



Features

- Neutral pH
- Alcohol Free
- Polymeric Biguanide (2000ppm)
- Non-ionic Surfactant
- 5th Generation Quaternary Ammonium (3000ppm)
- Dye Free
- Fragrance Free
- Non-corrosive
- Non-flammable



Atomic Ready is recommended for fixed surfaces and non-critical items:

Oximeters, BP cuffs, Stethoscopes, Thermometers, Oil cans, Bedpans, Buckets, Basins, Diuresis bottles, among others.



Oximeter



Stethoscopes



Blood Pressure Device



Thermometer



Bedpan



Oiler

How to use:

To clean and disinfect fixed surfaces, apply Atomic Ready directly to the area using a spray bottle, sprayer or a clean cloth dampened with the product.

To clean and disinfect non-critical items, wash with Atomic Ready, leave to act for 2 minutes, if necessary, rinse and dry the items.

When handling the product, we recommend the use of PPE.

